



## F1-16-MKII Impact of a Jet

Water is discharged vertically through a nozzle to strike a target carried on a stem, which extends through the cover. The dead weight of the moving parts are counterbalanced by a compression spring.

**PRODUCT CODE: F1-16-MKII**

**Categories:** Agricultural, Educational, F1-16-MKII Impact Of A Jet

**Tags:** Armfield , Compression Spring , Counterbalanced , Dead Weight , F1-16-MKII , Hemispherical Target , Impact Of A Jet , Momentum Theory , Nozzle , Vertical Force , Weights

### DESCRIPTION

Water is discharged vertically through a nozzle to strike a target carried on a stem, which extends through the cover. The dead weight of the moving parts are counterbalanced by a compression spring. The vertical force exerted on the target plate is measured by adding the weights supplied to the weight pan.

The apparatus consists of a cylindrical clear acrylic fabrication with provision for levelling. Water is fed through a nozzle and discharged vertically to strike a target carried on a stem, which extends through the cover. A weight carrier is mounted on the upper end of the stem.

The dead weight of the moving parts is counterbalanced by a compression spring. The vertical force exerted on the target plate is measured by adding the weights supplied to the weight pan until the mark on the weight pan corresponds with the level gauge.

A total of eight targets are provided.

### TECHNICAL SPECIFICATIONS

- Nozzle diameter: 8mm
- Distance between nozzle & target plate: 40mm
- Diameter of target plate: 36mm
- Target plates:
  - – 180° hemispherical target
  - – 120° target (cone)
  - – Flat target
  - – 30° target
  - – 60° target

- – Cup 135°
- – Oblique 30/150°
- – Oblique 45/135°

Requires Hydraulics Bench Service unit F1-10


## FEATURES & BENEFITS

### Demonstration Capabilities


- Principle of linear momentum
- To investigate the reaction forces produced by the change in momentum of a fluid flow
- Measurement of the forces produced by a jet impinging on solid surfaces which produce different degrees of flow deflection

### Practical Training Exercises

- Principle of linear momentum
- To investigate the reaction forces produced by the change in momentum of a fluid flow
- Measurement of the forces produced by a jet impinging on solid surfaces which produce different degrees of flow deflection



F1-16-MKII Impact of a Jet Graph



F1-16-MKII Impact of a Jet Diagram

## Dimensions

**Length:** 0.325m

**Width:** 0.20m

**Height:** 0.50m

## Order Codes

F1-16-MKII

## YOU MAY ALSO LIKE...



**F1-10  
HYDRAULICS  
BENCH**

[Read more](#)



**F1-33 PITOT  
TUBE  
DEMONSTRATO  
R**

[Read more](#)

.....  
Armfield Limited. Unit 10 Headlands Business Park, Salisbury Road, Ringwood, Hampshire BH24 3PB.  
England. Tel: +44(0) 1425 478781

