




F1-21-MKII Flow Meter Demonstration

This accessory is designed to introduce students to three basic types of flow meter:

- Venturi meter
 - Variable area flowmeter (Rotameter)
 - Orifice plate
- 8 pressure tapings are connected and displayed on the manometer bank to visualise pressure profiles

PRODUCT CODE: F1-21-MKII

Categories: Educational, F1 Series Fluid Mechanics

Tags: Armfield , F1-21 , F1-21-MKII , Flow Control Valve , Flow Measurement , Flow Meter Demonstration , Fluid Mechanics , Head Loss Characteristics , Hydraulics Bench , Manometers , Orifice Plate , Pressure Drops , Pressure Tappings , Variable-Area Flow Meter , Venturi Meter , Volumetric Measuring

DESCRIPTION

The equipment consists of a Venturi meter, variable area meter and orifice plate, installed in a series configuration to permit direct comparison. A flow control valve permits variation of the flow rate through the circuit. Pressure tapings are incorporated so that the head loss characteristics of each flow meter may be measured. These tapings are connected to an eight-tube manometer bank incorporating a manifold with an air bleed valve.

Pressurisation of the manometers is facilitated by a hand pump. The circuit and manometer are attached to a support framework, which stands on the working top of the Hydraulics Bench. The bench is used as the source of water supply and for volumetrically calibrating each flow meter.

- Optional Pitot Static Tube
- Optional Flow Nozzle

TECHNICAL SPECIFICATIONS

Manometer range: 0-400mm

Number of manometer: tubes 8

Orifice plate diameter: 17mm

Variable area meter: 2-20 l/min

VENTURI DIMENSIONS

Throat diameter: 14mm

Upstream pipe diameter: 26mm

Upstream taper: 21° inclusive
Downstream taper: 9° inclusive

Requires Hydraulics Bench Service unit F1-10/F1-10-2

FEATURES & BENEFITS

This accessory is designed to introduce students to three basic types of flow meter:

- Venturi meter
- Variable area flowmeter (Rotameter)
- Orifice plate
- 8 pressure tapings are connected and displayed on the manometer bank to visualise pressure profiles
- Optional Pitot Static Tube
- Optional Flow Nozzle

Experimental Content

- To investigate the operation and characteristics of a Venturi meter, variable area flowmeter and orifice plate including accuracy and energy losses
- Comparison of pressure drops across each flow measurement device
- Calibrating each flow meter using the volumetric measuring tank of the hydraulics bench
- Application of the Bernoulli equation for incompressible fluids
- Optionally available Pitot Static tube and Flow Nozzle

Video (English and Spanish)

Watch a short video about the F1-21-MKII in English:

Mira un breve video sobre el F1-21-MKII en español:



Requirements

F1-10 Hydraulics Bench to provide water flow.

Dimensions

Length: 0.68m

Width: 0.33m

Height: 0.83m

Order Codes

F1-21-MKII

YOU MAY ALSO LIKE...



**F1-15-MKII
BERNOULLI'S
THEOREM
DEMONSTRATI
ON**

[Read more](#)



**F1-16-MKII
IMPACT OF A
JET**

[Read more](#)



**C6-MKII-10 –
FLUID FRICTION
MEASUREMENT
S**

[Read more](#)



**F1-10
HYDRAULICS
BENCH**

[Read more](#)



**F1-22-MKII
ENERGY LOSSES
IN BENDS AND
FITTINGS**

[Read more](#)

Armfield Limited. Unit 10 Headlands Business Park, Salisbury Road, Ringwood, Hampshire BH24 3PB.
England. Tel: +44(0) 1425 478781