



FEL6 – Demonstration Lysimeter

Categories: Agricultural, Discipline / Sectors, Educational, Environmental, FE Series - Irrigation Water Management, FEL6 - Demonstration Lysimeter, Geology & Geography

Tags: Drainage , FE Series , Irrigation Water Management , Manometer , Precipitation

DESCRIPTION

The Armfield FEL6 lysimeter consists of three 300mm diameter containers in which a variety of crop types can be grown in any soil type.

Each container can then in turn be placed on a hydraulically mounted plate, which is used to monitor system weight changes arising from evapotranspiration, precipitation and drainage.

Despite its simplicity, the apparatus has a high degree of sensitivity and may be used for student project work as well as for demonstration purposes.

TECHNICAL SPECIFICATIONS

FEATURES & BENEFITS

Experimental content

- Study of the evapotranspiration measurement by use of a lysimeter
- Study of a lysimeter operation
- Determine lysimeter sensitivity to gains and losses in weight at different manometer angles and to relate this to gains or losses in water
- Demonstrate the components of the soil-plant-water system using the lysimeter
- Measuring the variations of water loss or gain from the soil and plant due to evapotranspiration, precipitation and drainage processes

- Demonstrate the principles of relationships between reference crop evapotranspiration, maximum and actual evapotranspiration
- Determination of the plants water usage
- Determine changes in diurnal evapotranspiration as measured by the lysimeter

Video

Ordering specification

Demonstration Lysimeter Apparatus comprising:

- Base and inner disc to support soil filled container and plant, resting on a hydraulic sensing device located in the base that is connected to a graduated water column mounted above the lysimeter
- Complete with spare plant containers
- Complete with instruction manual



Requirements

Cold water supply

Shipping Specifications

PACKED AND CRATED SHIPPING SPECIFICATIONS

Volume: 0.6m³

Gross Weight: 100kg

Dimensions

Length: 0.43m

Width: 0.312m

Height: 0.080m

Order Codes

FEL6

Armfield Limited. Unit 10 Headlands Business Park, Salisbury Road, Ringwood, Hampshire BH24 3PB.
England. Tel: +44(0) 1425 478781

