



FT110DP Benchtop Rapid Extractor



The Armfield FT110DP Benchtop Rapid Extractor is used to extract the active constituents from small (<1 litre) samples of materials.

This single piston pressure system gives rapid results from small volumes while giving flexibility to the operator with the processing parameters.

The three phases of operation are pressurisation, relaxation and dynamic percolation.

PRODUCT CODE: FT110DP

Categories: Armfield Edible And Essential Oils, Armfield Educational Food Technology, Armfield Flavours And Fragrances, Armfield Nutraceuticals, Armfield Pharmaceuticals, Chemical, FT Series - Educational Food Technology, FT110DP Rapid Extractor, IFT Industrial Food Technology

Tags: Agitation , Armfield , Dynamic Extraction , Extractor , FT110 , FT110DP , Ift , Industrial-Food-Tech , Percolation , Pressure , Rapid Extractor , Solvent Vessel , Speed-Machine , Static Extraction

DESCRIPTION

The Armfield FT110DP uses high pressure and a combination of both static and dynamic extraction phases to achieve a rapid extraction of the active materials with minimum degradation to the product. In the dynamic phases, the solvent is passed through the material providing a forced percolation and agitation.

The process pressure, number of cycles and length of each phase are programmable by the user. In this way the effect of the extraction parameters on the end product can be explored and an optimum set of parameters determined for each product.

The FT110DP vessel is loaded with solvent and primed manually. A single pneumatic piston is used to generate the pressure in the vessel and to provide the solvent percolation using a novel system of valves.

Included are two 5 litre external vessels. One for recirculation and solvent storage and one for solvent recovery.

ADDITIONAL INFORMATION

Voltage

USA, UK, Asia

TECHNICAL SPECIFICATIONS

Process vessel volume: 0.5L, 1L, 2L

Total solvent capacity: 5L

Maximum pressure: 8 bar

External Vessel Capacity: 5L each (two supplied as standard)

Number of Pistons: 2 Pistons

FEATURES & BENEFITS

- Fast response. Results in hours rather than days
- Small scale, only one litre of raw material is required
- Simple cleaning and maintenance procedures
- Uses ethanol, glycerol, water or similar safe polar solvents
- Reproducible results
- Fast, efficient solid/liquid extraction of active constituents from plants, herbs, fruits and other similar materials
- Computer Controlled
- Data Logging
- Accurate scale up
- Remote Control via VPN
- Touch Screen Membrane Control Panel
- High pressure, room temperature extraction, excellent for sensitive material
- Unaltered natural properties and characteristics of the active constituents, no discolouration or oxidation
- Easy to use: microprocessor control and real time display of sequencing and process parameters
- Automatic pressure control
- Automatic alarm indication of process anomalies
- Highly flexible extraction parameters
- Includes one process vessel, two solvent vessels

Applications

- Natural Food Flavours
- Fruit Flavours
- Herbs and Plant extracts
- Spirit liquors
- Tinctures and infusions
- Herbal extracts and homeopathy
- Cosmetics
- Pharmaceuticals



Requirements

- Electrical supply: Single phase (see ordering codes)
- Compressed Air (or N₂) 8bar

Shipping Specifications

PACKED AND CRATED SHIPPING SPECIFICATIONS

Volume: 0.3m³

Gross weight: 40Kg

Dimensions

Length: 0.70m

Width: 0.35m

Height: 0.55m

Order Codes

- **FT110DP-A:** 220-240V/1ph/50Hz, 0.5 kW
- **FT110DP-B:** 115V/1ph/60Hz, 0.5 kW
- **FT110DP-G:** 220-240V/1ph/60Hz, 0.5 KW

YOU MAY ALSO LIKE...



**FT18-MKII
MODULAR
CROSS FLOW
FILTRATION**

[Read more](#)



**FT25 SCRAPPED
SURFACE HEAT
EXCHANGER
SERIES**

[Read more](#)



**FT140X –
MIXING
VESSELS**

[Read more](#)

